

How to Disinfect and Use the Microchip Dispenser



Disinfect

1. Remove lid.
2. Leaving the round chip dispenser IN the bottle, pour a cold-sterile solution of your choice (alcohol, Nolvosan, etc.) into the bottle filling it to cover the dispenser.
3. Soak for 15 to 20 minutes. Do NOT leave the chips soaking indefinitely. You do not need to soak them every time you use the dispenser; just the first time.
4. Pour off liquid. Remove the dispenser from the bottle, allowing any excess fluid to drain out of the pin holes in the bottom of the dispenser.
5. Gently shake the chips to make sure all liquid is drained. Save the storage bottle.

Load the Chip

1. With the chip dispenser UPRIGHT, hold it firmly in one hand and rotate the top with your other hand until the loading-hole is lined up with a chip below. **NOTE:** Microchip damage can occur if you try to rotate the top while a chip is protruding from the cylinder.
2. Pull back the plunger rod of a Microchip Syringe (needle attached).
3. Slowly insert the needle into the open hole of the dispenser top. Holding the cylinder on a comfortable angle, you will now twist-slide the needle over the exposed chip. Snug the chip up into the needle. The dispenser is designed so that only one chip can be removed at a time, so don't worry about spillage.
4. Once the chip is loaded, remove the syringe from the dispenser. The chip should be positioned in the back of the needle, with only the tip of the chip visible at the opening of the bevel.
5. Be careful at this point not to prematurely bump or push on the plunger rod, as this will eject the chip.

Proceed with microchip injection following instructions provided with your products.

If you have any questions or concerns, please call (800) 434-2843

Microchip ID Systems, Inc.

www.microchipidsystems.com